


AIR-PUMP UNIT

The air puff unit is designed in compliance with,

- *Medical device **standards**
- *Ensuring it meets regulatory requirements for **safety**
- *Performance and quality control
- *The Air-puff **doesn't harm** the eye.

The air puff unit consists of

1. AIR-PUMP
2. PRESSURE REGULATOR
3. SOLENOID VALVE
4. NOZZLE
5. PIPE STRUCTURE



1. AIR-PUMP:

To generate the required puff of air ,a **12V DC vacuum air pump** will be used.

- *Rated current = **300mA**
- *Air pump rate = **4LPM**
- * Maximum pressure = **750mmHg**
- *Operation will be controlled by a **push button** and the **ESP32**
- *It has a length of 68mm and a diameter of 24mm,small in size



2. PRESSURE REGULATOR:

It regulates the pressure coming **from the air pump** and delivers it **to solenoid valve**.

- ***Miniature air pressure** regulator series
- *Ensures air puff is strong enough to flatten the cornea
- *Maximum supply pressure = **17.2 bar** (12901.06mmHg)
- *Can regulate up to desired air pressure (**50-55 mmHg**)
- *Regulate up to **0-4.1bar**



3. SOLENOID VALVE:

The solenoid valve **opens and closes within milliseconds**, providing the **quick air pulse** necessary for accurate tonometry.

The solenoid valve controls the release of a **precise and rapid air puff** directed toward the eye for IOP measurement.

- *Flow model = **unidirectional**
- *Response time-less than **50ms**
- *Operating pressure = **1.2Mpa**
- *Flow rate = **0.23Cv(fluid flow)**
- *Valve **type-2 way** normally open



4. NOZZLE:

It focuses air puff towards Center of the cornea (3-4mm), ensures the air is delivered in precise.

- *Customized 3D printed **ROUND-JET** model
- *Diameter = 2mm
- *Length = 2cm



From this Nozzle ,we can ensure that delivered **air-pressure** will be **55mmHg**.

5. PIPE STRUCTURE:

Ensuring a **smooth flow of air** within the system ,without **significant loss of pressure**.

Designed to maintain consistent pressure and velocity profiles for optimal performance in the air puff mechanism.

- *The pipe structure follows **BERNOULLI'S PRINCIPLE TO AVOID PRESSURE LOSS**.
- *Material used in the construction of pipe is **Polyurethane (PU)**.

Finally by **calibrating** all the units output ,we can get **controlled output (55mmhg)** to measure the **IOP accurately** without harming the user.

For orientation of cornea : Neon glittering paint is coated in the outer rim of nozzle , to achieve the correct orientation, X shaped mark is made in the hollow part of nozzle, so that the IR emitter and photodiode which are present in the sides of nozzle can align correctly.