•

ELECTRICAL UNIT

- The electrical unit plays a crucial role in controlling and managing the device's components, ensuring accurate and reliable intraocular pressure (IOP) measurement.
- It consists of
 - 1. Power supply
 - 2.Controller
 - 3. Display
 - 4.Transistor
 - 5.Mosfet
 - 6.Buttons

POWER SUPPLY:

Rechargeable lithium-ion battery

- Output voltage 12v
- Capacity -1200mAH
- Time taken to charge 24 minutes.
- Battery life 1 hour 30 minutes
- **Dimensions** 6.8 x 5.5 x 2 cm
- **Weight** 150g



CONTROLLER:

ESP32 WROOM

- Type ESP-32 with Wi-Fi 802.11 b/g/n
- Operating voltage 5V.
- Output voltage 3.3-5V.
- **Dimensions** 5 x 2.8 x 1.5 cm
- Weight 80g



ESP32 will be responsible for controlling the logic for the air pump, solenoid valve, IR emitter, and reading data from the photodiode.


DISPLAY:

- **Type** OLED Module
- Protocol 12C
- Support Voltage: 3.3 5V.
- Viewing angle >160°
- **Dimensions** 27mm x 27mm x 4mm
- Weight 10g

Delivers a better view of output to the user.

TRANSISTOR:

- It is used for switching and amplification purpose
- Switching: Control solenoid valves for air release
- Amplification: Boost weak sensor signals for better processing
- Control: Enable precise regulation of components like air pumps
- Using TIP102 transistor





5.MOSFET:

 The purpose of the MOSFET in this design is mainly related to power management, switching, and precise control of the electrical components such as the solenoid valve, air pump, or sensor circuitry.



Using MP1584 MOSFET

6.BUTTONS:

It consists of Power button and Push button

Medical-Grade Plastics:

ABS (Acrylonitrile Butadiene Styrene):

Properties: Durable, impact-resistant, and relatively lightweight.

Advantages: Easily molded into different shapes, cost-effective, and used in medical devices.

Applications: Commonly used in enclosures for medical instruments due to its strength and ease of processing.

Power button:

- It is used for on /off purpose to perform the operation
- Operating voltage:12V DC
- Operating current:2A MAX
- It is of **12mm**
- It is used to on the device and for releasing of air puff

Push button:

- It is used for adjusting the value of CCT
- It uses 2 push buttons of 5mm



